



External Implant System

*Value & Simplicity with
uncompromised Science*

BIOHORIZONS[®]
SCIENCE • INNOVATION • SERVICE



BIOHORIZONS[®]

external

Thread design demonstrates crestal bone preservation and 99.5% success rate in a 5-year prospective multi-center study.

Kline et al Implant Dent 2002;11:224-234

Higher bone-to-implant contact and reverse-torque than other thread designs.

Steigenga et al J Periodontol 2004;75:1233-1241

Predictable success in multiple Immediate Load studies.

Degidi and Piattelli J Periodontol 2003;74:225-241

Misch and Degidi Clin Implant Dent Relat Res 2003;5:17-28

Degidi et al Clin Implant Dent Relat Res 2006;8:123-134

3inOne : 3 functions / 1 component



Surgical Mount



Impression Coping
(with Ball-top Screw)



Esthetic Abutment

BioHorizons External provides a 3inOne[™] abutment with every implant.

What does your implant company give you?

- Biomechanical thread design maintains crestal bone levels and increases bone-to-implant contact.
- Esthetic gold-hued 3inOne Abutment lowers overhead and helps increase productivity.
- Implants and prosthetics are color-code matched for failsafe component mating.
- Surgical depth drills with built-in stops for added safety in critical areas.

Product Specifications

Material:	Titanium Alloy – Ti-6Al-4V
Surfaces:	Resorbable Blast Texturing (RBT) or Hydroxylapatite (HA)
Diameters:	Ø3.5mm Ø4.0mm Ø5.0mm Ø6.0mm
Lengths:	9mm 10.5mm 12mm 15mm
Prosthetic Options:	Comprehensive line of stock & custom-castable abutments
Warranty:	Lifetime unconditional

BioHorizons • 2300 Riverchase Center • Birmingham • Alabama 35244



For more information on BioHorizons External System,
contact the BioHorizons Customer Care department at
888-246-8338
or visit us on the web at
www.biohorizons.com

BioHorizons USA
888-246-8338

BioHorizons Canada
866-468-8338

BioHorizons Spain
+34 91 713 10 84

BioHorizons UK
+44 8700 620 550

BioHorizons Germany
+49 7661-909989-0

BioHorizons Australia
+61 2 8399 1520

BioHorizons Mexico
+52 55 5545 1297

BioHorizons Chile
+56 2 475 7230